

LPR 2000

LUMBAR PAIN RELIEVER

The **LPR 2000** deploys unique **opposing force** technology to accelerate healing in the lumbar spine. Patients often experience pain relief from the compression and postural improvement that the **LPR 2000** provides.



**Substantial Pain Relief
with Unparalleled
Comfort and Support**

- Relieves load on lumbar discs
- Disc space creation can alleviate nerve root impingement
- Low profile design for maximum patient compliance and comfort
- Unique design lifts and holds the patients' anterior weight
- Easy to clean and dry
- Posterior panel extends from the Sacrococcygeal Junction to the T-9 vertebrae
- Side lateral panels provide maximum lateral stability

American Medical Products LLC

LPR 2000

LUMBAR PAIN RELIEVER

SIZING

Measure around waist at navel

Size	Fits Waist
Small	23 - 27"
Medium	28 - 32"
Large	33 - 38"
X-Large	39 - 43"
XX-Large	44 - 48"
XXX-Large	49 - 54"
XXXX-Large	55 - 59"
XXXXX-Large	60 - 64"



How Does the LPR 2000 Work?

When the straps are pulled properly (outward and upward at a 45 degree angle), the straps work simultaneously to lift and compress the patients' anterior weight. The result is a distinct postural improvement and often a significant improvement in pain symptoms.

Rigid Posterior Panel

The interior rigid posterior panel is fabric covered and comfortably stabilizes the lumbar sacral region.

Indicated for:

- Spondylosis
- Chronic Lumbar Instability
- Post-Operative Rehabilitation
- Degenerative Disc Disease
- Chronic Back Pain
- Spinal Stenosis
- Herniated Lumbar Disc

LSO, Sagittal-Coronal control, with rigid anterior and posterior panels, posterior extends from sacrococcygeal junction to T-9 vertebrae, lateral strength provided by rigid lateral frame/panels, produces intracavitary pressure to reduce load on intervertebral discs, includes straps, pendulous abdomen design, prefabricated, includes fitting and adjustment.

American Medical Products LLC

151 B Industrial Way East • Eatontown, NJ 07724
Tel: 888.487.3101 • Fax: 888.487.7016
Email: info@americanmedicalproducts.com